


Gastight Syringes | 1800 Series



The 1800 series syringe is the Gastight version of the 800 series syringes. The 1800 series was designed to eliminate the possibility of plunger damage. It has the same liquid handling capabilities as the 1700 series but the extended plunger eliminates breakage and allows the user some control over the dispense speed.

Recommended Use: This syringe is intended for applications where the plunger on a 1700 series syringe might be bent.

Accessories: Reproducibility Adapter
Maximum Temperature: Cemented Needle: 122° F (50 ° C)
 Removable Needle: 239° F (115° C)
Maximum Pressure: 1,000 PSI, except 250 µL: 500 PSI

1800 SERIES ORDERING INFORMATION								
		VOLUME	10 µL	25 µL	50 µL	100 µL	250 µL	
		MODEL	1801	1802	1805	1810	1825	
		GAUGE	26s	22s	22s	22s	22s	
Removable Needle Syringes & Repl Needles		Removable Needle (RN) Syringe (point style 2)	84877	84880	84883	84886	84889	
		Removable Needle (RN) Syringe (needle sold separately)	7659-01	7660-01	7661-01	7662-01	7663-01	
Removable Needle Syringes & Repl Needles		Replacement Needles (sold separately)	Point Style	Small RN	Small RN	Small RN	Small RN	Large RN
			2	7758-02	7758-03	7758-03	7758-03	7779-03
			3	7768-01	7770-01	7770-01	7770-01	7780-03
			5	7784-07	7784-05	7784-05	7784-05	7784-01
ACCESSORIES (Syringe not included)								
Reproducibility (Chaney) Adapter			32146	32146	32146	32146	32146	

For Syringe Calibration Service, add the prefix "CAL" to the beginning of the syringe's part number.

Note: Needles are 51 mm

* "s" indicates a thicker wall