

Microliter Syringes | 900 Series



The 900 series has the same liquid handling capabilities as the 800 series but is a more economical version. The 900 series does not have any replaceable parts making the initial cost less than the 800 series.

Recommended Use: This syringe is perfect when a robust syringe is desired and budget is a consideration.
Accessories: Reproducibility Adapter

Accessories: Reproducibility Adapter
Maximum Temperature: Cemented Needle: 122° F (50 ° C)
 Removable Needle: 239° F (115° C)
Maximum Pressure: 2,000 PSI

900 SERIES ORDERING INFORMATION

		VOLUME	5 µL	10 µL															
		MODEL	95	901															
		GAUGE	26s	26s															
		Cemented Needle (N) Syringe (point style 2)	87920	80360															
Removable Needle Syringes & Repl. Needles		Removable Needle (RN) Syringe (point style 2)	87925	80370															
		Removable Needle (RN) Syringe (needle sold separately)	7647-01	7648-01															
		<table border="1"> <thead> <tr> <th>Replacement Needles (sold separately)</th> <th>Point Style</th> <th>Small RN</th> <th>Small RN</th> </tr> </thead> <tbody> <tr> <td></td> <td>2</td> <td>7758-02</td> <td>7758-02</td> </tr> <tr> <td></td> <td>3</td> <td>7768-01</td> <td>7768-01</td> </tr> <tr> <td></td> <td>5</td> <td>7784-07</td> <td>7784-07</td> </tr> </tbody> </table>	Replacement Needles (sold separately)	Point Style	Small RN	Small RN		2	7758-02	7758-02		3	7768-01	7768-01		5	7784-07	7784-07	
Replacement Needles (sold separately)	Point Style	Small RN	Small RN																
	2	7758-02	7758-02																
	3	7768-01	7768-01																
	5	7784-07	7784-07																
ACCESSORIES (Syringe not included)																			
Reproducibility (Chaney) Adapter			32148	32148															

For Syringe Calibration Service, add the prefix "CAL" to the beginning of the syringe's part number.
 Note: Needles are 51 mm