

Do you know?

Sonic Reservoir Sensor

SR-13



Safe and Continuous Monitoring of Reagent Depletion and Waste Overflow

Notify waste level before overflow

Designed for safe handling of waste organic solvents

Monitor reagent and waste levels in real time!


Prevent LC system failure when solvent levels become critical

SR-13 Specifications

Sensing mechanism	Sonic
Type of sensor	2 types, for reagent and waste reservoirs
Number of sensors	Max. 3 ch. (selectable sensor type)
Alert function	LED (real time monitoring), Error buzzer and signal output (0-5 v)

FLOM Coporation

5-32-10 Shinmachi, Ome, Tokyo 198-0024 Japan
 Phone: +81-428-30-7454
 Website: www.flom.co.jp
 Email: info@flom.co.jp



INACOM
INSTRUMENTS

Tel. 0318-521151*
 E-mail: info@inacom.nl

any